
# Tip StackPack Reload System

Contamination-free handling, compact and reusable



## Contamination-free stacking system

Tips are preloaded in reusable, stacked 96-tip trays. The trays are designed to safely cover each individual 96-tip layer, shielding the preloaded tips in the stack like a protective screen.

- Certified PCR Performance Tested: Free of human DNA, DNase, RNase and PCR inhibitors
- Shielding function of stacked trays ensures contamination-free handling of preloaded tips

## Considerably reduced waste volume compared to tip boxes

The lid, trays and base section of the StackPack system enable multiple use. As a result, composite materials and storage space are reduced by more than 50%. All materials are fully recyclable.

- Easy transfer of individual layers of 96-tip trays into the Tip SystemBox workstation
- · Safety labels on both sides can be easily perforated to access the next layer of tips
- Tips can also be used directly from the stack

Order No.	Description	Volume	Colour code	Option	Packaging	
70.760.452		200 µl			480/StackPack 1,920/case	
70.760.451		250 µl			480/StackPack 1,920/case	
70.1130.460		10 µl			576/StackPack 2,304/case	

#### Tip StackPack in PCR Performance Tested quality

#### Tip SystemBox

Order No.	Description	Packaging/Case
95.1298.001	Box with hinged lid (w/o pipette tips), autoclavable*	6 pcs.

SARSTEDT AG & Co. P.O. Box 12 20 D-51582 Nümbrecht Phone +49 2293 3050 Fax +49 2293 305-122 info@sarstedt.com www.sarstedt.com \* Products made from PP can be autoclaved up to 121°C without significantly impairing their mechanical properties. Users should test for themselves if autoclaving has any other effect on individual applications.



NEW